

EXTERNAL MINI

	Hexagon	2.41
	Lower External Dia	Φ 3.5
	Body External Dia	Φ 4.1
	Color	2
	Depth of Hexagon	0.8
	Length	10

compatible : Nobel, US, Megagen, Dio, Dentis, Protect

	Hexagon	2.41
	Lower External Dia	Φ 3.5
	Body External Dia	Φ 4.1
	Color	2
	Depth of Hexagon	0.8
	Length	8.5

compatible : Nobel, US, Megagen, Dio, Dentis, Protect

	Hexagon	2.41
	Lower External Dia	Φ 3.3
	Body External Dia	Φ 5
	Color	5
	Depth of Hexagon	1.1
	Length	10.5

compatible : Seojeon(Life Core)

	Hexagon	2.41
	Lower External Dia	Φ 3.3
	Body External Dia	Φ 5
	Color	5
	Depth of Hexagon	1.1
	Length	10.5

compatible : Nobel, US, Megagen, Dio, Dentis, Protect

EXTERNAL WIDE

	Hexagon	3.41
	Lower External Dia	Φ 5.1
	Body External Dia	Φ 6
	Color	3
	Depth of Hexagon	1.1
	Length	8.5

compatible : Nobel, US, Endopore, Dio, Dentis, Protect, oneplant

	Hexagon	3.41
	Lower External Dia	Φ 5.1
	Body External Dia	Φ 6
	Color	5
	Depth of Hexagon	0.8
	Length	10.5

compatible : Nobel, US, Endopore, Dio, Dentis, Protect

	Hexagon	3.41
	Lower External Dia	Φ 5.1
	Body External Dia	Φ 7
	Color	3
	Depth of Hexagon	0.8
	Length	8.5

compatible : Nobel, US, Endopore, Dio, Dentis, Protect

	Hexagon	3.41
	Lower External Dia	Φ 5.1
	Body External Dia	Φ 7
	Color	3
	Depth of Hexagon	0.8
	Length	10.5

compatible : Nobel, US, Endopore, Dio, Dentis, Protect

	Hexagon	3.41
	Lower External Dia	Φ 5.0
	Body External Dia	Φ 6
	Color	4
	Depth of Hexagon	0.8
	Length	12

compatible : Nobel, US, Endopore, Dio, Dentis, Protect

	Hexagon	3.41
	Lower External Dia	Φ 5.0
	Body External Dia	Φ 6
	Color	4
	Depth of Hexagon	0.8
	Length	12

compatible : 3R, Megagen, Cowellmedi

	Hexagon	3.41
	Lower External Dia	Φ 5.0
	Body External Dia	Φ 6
	Color	4
	Depth of Hexagon	1.1
	Length	12

compatible : Biolock

	Hexagon	3.41
	Lower External Dia	Φ 5.0
	Body External Dia	Φ 6
	Color	4
	Depth of Hexagon	1.1
	Length	12

compatible : Biolock

EXTERNAL REGULAR

	Hexagon	2.71
	Lower External Dia	Φ 4.1
	Upper Angle	15
	Color	-
	Depth of Hexagon	0.8
	Length	9.74

compatible : Nobel, 3I, Seojeon, Endopore, US, Megagen, Dio, Dentis, Protect, Cowellmedi

	Hexagon	2.71
	Lower External Dia	Φ 4.1
	Upper Angle	25
	Color	-
	Depth of Hexagon	0.8
	Length	9.42

compatible : Nobel, 3I, Seojeon, Endopore, US, Megagen, Dio, Dentis, Protect, Cowellmedi

	Hexagon	2.71
	Lower External Dia	Φ 4.1
	Body External Dia	Φ 5
	Color	3
	Depth of Hexagon	0.8
	Length	8.5

compatible : Nobel, 3I, Seojeon, Endopore, US, Megagen, Dio, Dentis, Protect, Cowellmedi

	Hexagon	2.71
	Lower External Dia	Φ 4.1
	Body External Dia	Φ 5
	Color	5
	Depth of Hexagon	0.8
	Length	10.5

compatible : Nobel, 3I, Seojeon, Endopore, US, Megagen, Dio, Dentis, Protect, Cowellmedi

	Hexagon	2.71
	Lower External Dia	Φ 4.1
	Body External Dia	Φ 6
	Color	3
	Depth of Hexagon	0.8
	Length	8.5

compatible : Nobel, 3I, Seojeon, Endopore, US, Megagen, Dio, Dentis, Protect, Cowellmedi

	Hexagon	2.71
	Lower External Dia	Φ 4.1
	Body External Dia	Φ 6
	Color	5
	Depth of Hexagon	0.8
	Length	10.5

compatible : Nobel, 3I, Seojeon, Endopore, US, Megagen, Dio, Dentis, Protect, Cowellmedi

	Hexagon	2.71
	Lower External Dia	Φ 4.1
	Body External Dia	Φ 6
	Color	4
	Depth of Hexagon	0.8
	Length	12

compatible : Nobel, 3I, Seojeon, Endopore, US, Megagen, Dio, Dentis, Protect, Cowellmedi

	Hexagon	2.71
	Lower External Dia	Φ 4.1
	Body External Dia	Φ 6
	Color	4
	Depth of Hexagon	0.8
	Length	12

compatible : Nobel, 3I, Seojeon, Endopore, US, Megagen, Dio, Dentis, Protect, Cowellmedi

	Hexagon	2.71
	Lower External Dia	Φ 4.1
	Body External Dia	Φ 9
	Color	4
	Depth of Hexagon	0.8
	Length	9.8

compatible : Nobel, 3I, Seojeon, Endopore, US, Megagen, Dio, Dentis, Protect, Cowellmedi

	Hexagon	2.71
	Lower External Dia	Φ 4.1
	Body External Dia	Φ 5
	Color	4
	Depth of Hexagon	1.1
	Length	8.5

compatible : oneplant, Biolock

	Hexagon	2.71
	Lower External Dia	Φ 4.1
	Body External Dia	Φ 5
	Color	5
	Depth of Hexagon	1.1
	Length	10.5

compatible : oneplant, Biolock

	Hexagon	2.71
	Lower External Dia	Φ 4.1
	Body External Dia	Φ 6
	Color	3
	Depth of Hexagon	1.1
	Length	8.5

compatible : oneplant, Biolock


	Hexagon	2.71
	Lower External Dia	Φ 4.1
	Body External Dia	Φ 6
	Color	5
	Depth of Hexagon	1.1
	Length	10.5

compatible : oneplant, Biolock


	Hexagon	2.71
	Lower External Dia	Φ 4.1
	Body External Dia	Φ 6
	Color	4
	Depth of Hexagon	1.1
	Length	12

compatible : oneplant, Biolock


INTERNAL ASTRA TYPE

	Lower External Dia	Φ 4.1
	Body External Dia.	Φ 4.5
	Color	3.3
	Length	13


compatible : Astra, Implantium, GS II, Megafix, Dentis

	Lower External Dia	Φ 4.1
	Upper Angle	Φ 5.5
	Color	3.3
	Length	13


compatible : Astra, Implantium, GS II, Megafix, Dentis

	Lower External Dia	Φ 4.1
	Body External Dia	Φ 5.5
	Color	4.8
	Length	14.5


compatible : Astra, Implantium, GS II, Megafix, Dentis

	Lower External Dia	Φ 4.1
	Body External Dia	Φ 6.5
	Color	3.3
	Length	13


compatible : Astra, Implantium, GS II, Megafix, Dentis

	Lower External Dia	Φ 4.1
	Body External Dia	Φ 6.5
	Color	4.8
	Length	14.5


compatible : Astra, Implantium, GS II, Megafix, Dentis

	Lower External Dia	Φ 4.1
	Body External Dia	Φ 7.5
	Color	4.8
	Length	14.5


compatible : Astra, Implantium, GS II, Megafix, Dentis

	Lower External Dia	Φ 4.1
	Body External Dia	Φ 4.5
	Upper External Dia.	Φ 4.2
	Length	13


compatible : Astra, Implantium, GS II, Megafix, Dentis

	Lower External Dia	Φ 4.5
	Body External Dia	Φ 5.5
	Color	3.3
	Length	13


compatible : Astra(4.5)

	Lower External Dia	Φ 4.1
	Body External Dia	Φ 5.5
	Upper External Dia.	Φ 4.9
	Length	14.5


compatible : Astra, Implantium, GS II, Megafix, Dentis

	Lower External Dia	Φ 4.1
	Body External Dia	Φ 5.5
	Color	4.8
	Length	14.5


compatible : Astra(4.5)

	Lower External Dia	Φ 4.1
	Body External Dia	Φ 6.5
	Upper External Dia.	Φ 4.9
	Length	14.5


compatible : Astra, Implantium, GS II, Megafix, Dentis

	Lower External Dia	Φ 4.1
	Body External Dia	Φ 7.5
	Upper External Dia.	Φ 5.7
	Length	14.5

compatible : Astra, Implantium, GS II, Megafix, Dentis


	Lower External Dia	Φ 4.1
	Body External Dia	Φ 5.5
	Upper Angle	17
	Length	14

compatible : Astra, Implantium, GS II, Megafix, Dentis


	Lower External Dia	Φ 4.1
	Body External Dia	Φ 5.5
	Upper Angle	17
	Length	12

compatible : Astra, Implantium, GS II, Megafix, Dentis


INTERNAL SUBMERGE TYPE

	Lower External Dia	Φ 4.1
	Body External Dia	Φ 6
	Color	4
	Length	11


compatible : 3I regular

	Lower External Dia	Φ 4.1
	Body External Dia	Φ 6
	Color	4
	Length	11


compatible : Zimmer (TSV)

	Lower External Dia	Φ 5.7
	Body External Dia	Φ 6
	Color	11
	Length	11


compatible : Zimmer (TSV)

	Lower External Dia	Φ 3.7
	Body External Dia	Φ 4
	Upper Angle	2
	Length	11


compatible : POIEX37

	Lower External Dia	Φ 4.1
	Body External Dia	Φ 5
	Color	1
	Length	10


compatible : 3I regular

	Lower External Dia	Φ 4.5
	Body External Dia	Φ 4.7
	Color	11
	Length	11

compatible : Zimmer (TSV)

	Lower External Dia	Φ 3.0
	Body External Dia	Φ 3.8
	Color	8
	Length	8

compatible : ITI

	Lower External Dia	Φ 3.7
	Body External Dia	Φ 4.5
	Upper External Dia.	Φ 2
	Length	11

compatible : POIEX42