

Gutta Percha Points



Made from the Finest Gutta Percha Uniformly Hand Rolled

- ✓ Combination of stiff and flexible properties provide ideal working qualities
- ✓ Uniform and non-distorting crimp
- ✓ Excellent radiopacity
- ✓ Long shelf life and freshness assured
- ✓ Suitable for vertical and lateral condensation
- ✓ Variety in the selection of sizes and packaging available






Special Length Marked Gutta Points



For Accurate Depth Measurement and Insertion (Patent Pending)

- ✓ Uniform and non-distorting crimp
- ✓ Variable sizes and packages available
- ✓ Non-Standardized (Accessory Sizes)
- ✓ Accurate insertion of the working length
- ✓ No need to mark working length on the points
(Assurance of the accurate insertion after filling)
- ✓ Time saving with safe, easy and fast filling
- ✓ Variety in the selection of sizes and packaging available

Available Packaging:

	<p>Spill proof (6 Compartment Dispenser)</p> <p>Packaging Code: P.21 Points: 120 points Order Code: GPCB L x W x H: 80x51x16(mm) Color Coded: Yes Length Marked: No Sizes: Standard & Non Standard</p>
	<p>Vials (6 Round Vials)</p> <p>Packaging Code: P.31 Points: 120 points Order Code: GPCV L x W x H: 80x51x16(mm) Color Coded: Yes Length Marked: No Sizes: Standard & Non Standard</p>
	<p>.04 Taper Gutta Percha</p> <p>Packaging Code: P.21 Points: 60 points Order Code: GP04 L x W x H: 80x51x16(mm) Color Coded: Yes Length Marked: Yes Sizes: ISO #15 - 80</p>
	<p>.06 Taper Gutta Percha</p> <p>Packaging Code: P.21 Points: 60 points Order Code: GP06 L x W x H: 80x51x16(mm) Color Coded: Yes Length Marked: Yes Sizes: ISO #15 - 80</p>
	<p>Endo Organizer (Spill proof 12 Compartments)</p> <p>Packaging Code: P.11 Points: 240 points Order Code: ENDOGP L x W x H: 110x110x30(mm) Color Coded: Yes Length Marked: No Sizes: ISO #15/40 (assorted only) ISO #15/80 (assorted only)</p>

Available Sizes:

ISO STANDARD

#8	●	#55	●
#10	●	#60	●
#15	○	#70	●
#20	●	#80	●
#25	●	#90	○
#30	●	#100	●
#35	●	#110	●
#40	●	#120	●
#45	○	#130	●
#50	●	#140	●

NON-STANDARD

- X-Fine
- Fine-Fine
- Medium-Fine
- Fine
- Fine-Medium
- Medium
- Medium-Large
- Large
- X-Large

ASSORTED

- 15/40
- 45/80
- Assorted (XF-M)