

Absorbent Paperpoints



Hand Rolled for Extra Absorbency

- ✓ Rigid enough to insert in canals without curling or bending
- ✓ Rolled firmly to hold shape when withdrawn from the canal
- ✓ Highly absorbent
- ✓ Sterilized individual packages
- ✓ Variety sizes and packages available

Available Sizes

ISO STANDARD

#8	-----	●	#55	-----	●
#10	-----	●	#60	-----	●
#15	-----	○	#70	-----	●
#20	-----	●	#80	-----	●
#25	-----	●	#90	-----	○
#30	-----	●	#100	-----	●
#35	-----	●	#110	-----	●
#40	-----	●	#120	-----	●
#45	-----	○	#130	-----	●
#50	-----	●	#140	-----	●








NON-STANDARD

- XX-Fine
- X-Fine
- Fine
- Medium
- Coarse
- X-Coarse

ASSORTED

- 15/40
- 45/80
- Assorted (XXF – XC)

Available Packaging

	<p>Spill proof (6 Compartment Dispenser)</p> <p>Packaging Code: P.21 Points: 200 points Order Code: PPCB L x W x H: 84x51x13(mm) Color Code: Yes Sizes: Standard & Non-Standard</p>
	<p>Vials (6 Round Vials)</p> <p>Packaging Code: P.31 Points: 200 points Order Code: PPCV L x W x H: 84x51x16(mm) Color Code: Yes Sizes: Standard & Non-Standard</p>
	<p>.04 Taper (Spill proof 6 Compartments)</p> <p>Packaging Code: P.21 Points: 60 Points Order Code: PP04 L x W x H: 84x51x16(mm) Color Code: Yes Sizes: ISO #15 - 80</p>
	<p>.06 Taper (Spill proof 6 Compartments)</p> <p>Packaging Code: P.21 Points: 60 Points Order Code: PP06 L x W x H: 84x51x16(mm) Color Code: Yes Sizes: ISO #15 - 80</p>
	<p>Cell Pack 180</p> <p>Packaging Code: P.61 Points: 180 points Order Code: PP180 L x W x H: 120x54x12(mm) Color Code: No Sizes: Non-Standard</p>
	<p>Cell Pack 200</p> <p>Packaging Code: P.61 Points: 200 points Order Code: PP200 L x W x H: 162x54x11(mm) Color Code: No Sizes: Standard & Non-Standard</p>
	<p>Endo Organizer (Spill proof 12 Compartments)</p> <p>Packaging Code: P.11 Points: 480 points Order Code: ENDOPP L x W x H: 110x110x30(mm) Color Code: Yes Sizes: ISO #15/40 (assorted only) ISO #15/80 (assorted only)</p>