



META DENTAL CORP.

BesQual S800 - Multi Wax Pot



- ✓ Individual temperature setting for each pot
- ✓ LCD Digital display for each pot
- ✓ Built-in memory function to keep previous temperature.
- ✓ Built-in thermostat accurately maintain desired temperature

ITEM CODE: BQ-S800
TEMPERATURE: 50°C ~ 120°C | 122°F ~ 248°F
WEIGHT: 0.69 Kg | 1.52 lb
DIMENSION (mm): 200 x 146 x 70